



Perfect Dent HENRY SCHEIN®

HENRY SCHEIN®

DENTAL DAM

• Natural Rubber Latex • Mint

ES Dique de goma • latex puro • menta

FR Digue dentaire • pur latex • menthe

DE HS-Kofferdam • natürlicher Latex • mint

IT Diga dentale • puro lattice di gomma • menta



NL Kofferdam • natuurlijk latex • mint

PT Dique dentário • látex puro de borracha natural • hortelã

EL Καουτσούκ φράγμα • φυσικό λατέξ • μέντα

TR Baraj • doğal lateks • nane

CS Dental Dam • přírodní latex • mátyový

HR Kofferdam • prirodni latex • mint

HU Kofferdam • természetes latex • menta

PL Koferdam stomatologiczny

RO Folie diga dentara • Latex

152x152mm (6" x 6")

36

POWDER 5°C / 41°F



HENRY SCHEIN

Rubber Dam & Clamps

HENRY SCHEIN

Rubber Dam

Natural Rubber Latex

Color: Green, Blue

Mint

152×152 mm (6"×6")

36 Sheets in pack

Thickness: Heavy (Providing great adaptation around teeth & does not tear easily)



Isolation
AS IT SHOULD BE

HENRY SCHEIN

Clamps

High-quality, stainless steel rubber dam clamps **PROVIDE:**

Nonreflective matte finish for better visibility

Individually labeled



Strong retention












Durable

Patient comfort

HENRY SCHEIN

Clamps

Henry Schein Clamps Set

Type	Description	Photo	
Winged	2	For bicuspid (Premolars) with flat beaks	
	2A	For bicuspid (Premolars) with flat beaks Similar to #2, but with larger beaks. For bicuspid	
	7	For big lower molars with flat beaks.	
	8A	For molars with deeply festooned beaks that can grab subgingival part	
	9	Most popular winged anterior clamp	
	12A	Especially for upper left and lower right molars with serrated beaks, paired with 13A	
	13A	Especially for lower left and upper right molars with serrated beaks, paired with 12A.	
	14	For partially erupted molars with deeply festooned beaks.	
	14A	For partially erupted molars with deeply festooned beaks. Larger than 14	
Wingless	W3	For small lower molars with flat beaks.	
	W8A	For molars with deeply festooned beaks that can grab subgingival part Similar to 8A but wingless.	

CONTACT US



perfectdenteg.com



perfectdenteg



perfectdenteg



info@perfectdenteg.com



+201552974571

+201018436971



74 Ammar ibn Yasser, above Metro Market,
Masr Al Jadidah, Cairo, Egypt