



DENTAL SENSOR

IMAGING INNOVATION



Powered
by



Perfect Dent



Size 1.3 full range
active area



Smart
Multi-functional
Imaging System



Waterproof IPX7



Swift and clear
imaging



Years
Warranty



Specification:



Durable bending times



Robust body



Ultra-slim and Ergonomic Design



Pixel Size: 18.5 μ m




Active Pixel Array (pixels): 1200 x 1600



Active Area: 29.6 x 22.2 mm





Imaging Time: 3 - 4 S



Outstanding Imaging Quality

CSI(Cesium Lodide) and CMOS technologies allows you to acquire high-definition images at lower radiation doses.



Contact Us



perfectdenteg.com



perfectdenteg



perfectdenteg



info@perfectdenteg.com



+201552974571
+201018436971



74 Ammar ibn Yasser, above Metro Market,
Masr Al Jadidah, Cairo, Egypt