

Powered
by



PERFECT
DENT



BLACK HAWK ENDOMOTOR


2
YEARS
WARRANTY

PD MAKING
THE FUTURE

RED IS TAKING
OVER EGYPT

SMART DESIGN

▪ **Contra:**

1. Mini head, slim neck
2. Reinforced Aluminum coated with corrosion resistant layer
3. Royal black color
4. **1:1** smooth contra provides the exact speed and torque of Motor
5.  rotatable contra angle



ERGONOMIC DESIGN

BLACK X BLACK PIECE OF ART



Weight
only 125g



Speed
150- 1000 rpm



Torque
0.6 - 5 N.cm



10 programs



MICRO USB type C,
rechargeable, fast charging



Powerful Battery
1200 mAh



Motions

1. FW (forward)
2. Rev (reverse)
3. OTR (optimum torque reverse)
4. REC (reciprocation)

*Fully adjustable angles

*Speed starts from 250-500 rpm

L/R

Suitable for
right & left handed operators




CONTACT US



perfectdenteg.com



perfectdenteg



perfectdenteg



info@perfectdenteg.com



+201552974571

+201018436971



74 Ammar ibn Yasser, above Metro Market,
Masr Al Jadidah, Cairo, Egypt