

# LAFOMED

## Autoclave

MAKE IT SAFER



#Autoclaves\_Masters

#أوتوكلافك\_بقي\_من\_بيرفكت

#التعقيم\_في\_أمان

WARRANTY  
2 YEARS

# LAFOMED



**Safe & Efficient**



**Built-in Drying**



**PERFECT  
DENT**

# LAFOMED




**23 Litre**



**5 Programs**



**PERFECT  
DENT**




**User-Friendly**




**Reliable Brand**



**PERFECT  
DENT**



|   |           |   |                  |    |                            |    |              |
|---|-----------|---|------------------|----|----------------------------|----|--------------|
| 1 | Screen    | 5 | Sterile chamber  | 9  | Outlet of clean water tank | 13 | Air switch   |
| 2 | Key       | 6 | Air filter       | 10 | Outlet of Waste water tank | 14 | Power switch |
| 3 | Handle    | 7 | Clean water tank | 11 | Waste drainage outlet      | 15 | SD card port |
| 4 | Door Seal | 8 | Safety valve     | 12 | Power cord                 | 16 | Printer      |



## Chamber Specification

**Dimensions (DN \* Length mm) :  $\Phi 254*450$**

## Parameter of Autoclave:

**Max working temperature: 137 °C**

**Max working pressure: 0.24 MPa**

**Safety valve set-pressure: 0.24 MPa**

**Safety valve open-pressure: 0.24~0.26 MPa**

**Capacity of clean water tank: 4 L**

**Capacity of waste water tank: 2.6 L**

**Power input: 220 V $\pm$ 22 V, 50 Hz, 2300 VA**

**Sound level:  $\leq 70$  dB(A)**

## Accessories:



**3 Instrument Trays**

**2.5 Kg**  
**Single layer**  
**maximum wight**



**Instrument**  
**Tray Handle**



**Drain Pipe**





# Perfect Dent



**Now, Best Autoclave has Advantage  
of your favourite company's customer service  
and extended warranty**





# Contact Us

---



perfectdenteg.com



perfectdenteg



perfectdenteg



info@perfectdenteg.com



+201552974571  
+201018436971



74 Ammar ibn Yasser, above Metro Market,  
Masr Al Jadidah, Cairo, Egypt