

**UDG**  
UNITED DENTAL GROUP








# M3 Pro Gold Files

- Root Canals Files
- Non cutting tip
- Sharp cutting edge assures its excellent resistance and flexibility at the same time
- **Certification:** CE, ISO

**M3-Pro Gold** is a series of multiple sequence files that applied with the most advanced metallurgy technology. The gold heat treatment enhances the flexibility and resistance to cyclic fatigue. Its **convex triangular cross-section** increase the cutting efficiency and preserve the detin as well. It can cover the needs of 95% clinical cases especially those curved canals.

## Non-cutting Tip

The **non-cutting** tip design allows each instrument to safely follow the secured portion of the canal while the small flat area on the tip enhances its ability to find its way through soft tissue and debris.



## High Flexibility

**M3-Pro Gold** proprietary processing allows it to follow canal curvature. As such, the file may retain the shape of the canal upon removal. This is normal. The file can be easily straightened or further curved.

## Polishing Technique

Advanced **triple** polishing technique, **Eliminate** stress point on the surface





# Gutta Percha Points

- Standard **Gutta Percha Points** are faster and effective
- Used for quick obturation after root canal preparation
- Excellent in workability with **non-distorted** crimp



# Sterile Absorbent Paper Points

- Long shelf life & Freshness Assured
- Stiff but flexible to provide ideal workability
- **Absorbent Paper Points** are faster and effective
- Used for quick obturation after root canal preparation





# CONTACT US

---



perfectdenteg.com



perfectdenteg



perfectdenteg



info@perfectdenteg.com



+201552974571

+201018436971



74 Ammar ibn Yasser, above Metro Market,  
Masr Al Jadidah, Cairo, Egypt